



Please contact Invacare EU Export Customer Services on :-
Tel : +33 247 626980 ; Email : serviceclient_export@invacare.com
Or visit our website : www.invacare.eu.com



alber[®] via GO[®]
Power add on for manual wheelchairs







Date :
Purchase order:
Customer code :
Distributor stamp :



- Descriptif :**
- New electric propulsion assistance: **via GO**
 - Very light : 5,85 kg
 - Remote to adjust speed and active the assistance
 - Compatible with 80% manual wheelchairs
 - Compatible with aircraft
 - Powerfull motor ,perform in slopes
 - Range : jusqu'à 16 km
 - Max user weight: 150 kg
 - Assembly possible par all dealers or **Invacare**

☐ Sélectionner les options choisies par une croix dans les cases ☒ Configuration standard ☐ Options sans supplément de prix

| | | | Euros | |
|---------------------|---|---|--------------|---|
| | | | Public price | |
| CODE | DESCRIPTION | NC : Not Compatible NC : 9999 : not compatible with reference 9999 | HT | |
| via GO V24 | | | | |
| 1595002 | <input type="checkbox"/> alber [®] via GO V24 - 5,5 km/h - <i>Charger and remote to adjust speed and ON/OFF .</i> | | 2 701,42 |  |
| OPTIONS IMPORTANTES | | | | |
| 1595020 | <input type="checkbox"/> Standard bracket for mounting on standard and lightweight wheelchairs with tube diameter 25 mm (can be mounted by all dealers) | | 236,97 |  |
| 1595019 | <input type="checkbox"/> Bracket set for mounting on wheelchairs with different tube diameter, comfort chairs and seat shell frame | | 284,36 |  |
| 1595080 | <input type="checkbox"/> Counterholder bracket control unit incl. screws for wheelchairs without drum brake | | 28,44 |  |
| 1584856-SM | <input type="checkbox"/> Installation of the complete fixing kit for Viago by Invacare (Not mandatory) (return shipping costs not included) Average assembly time: 5 working days upon receipt of the chair (In the case of an existing chair, plan to return it to our workshops within 10 days). | | 109,95 | |

INFORMATIONS CONCERNANT VOTRE FAUTEUIL :


Serial number: _____

Seat width (in cm): _____

Seat height (in cm):: _____

Distance between the ends of the quick-release axles (in cm, for folding wheelchair): _____

| | |
|--------------------------|-----------------------------------|
| <input type="checkbox"/> | Equipped with drum brakes |
| <input type="checkbox"/> | Equipped with antitippers |
| <input type="checkbox"/> | Equipped with transportation kit |
| <input type="checkbox"/> | Equipped with wheelbase extension |



Test the compatibility of the manual wheelchair using the DIP